Instruction manual

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WELLTEC OCM10

Oxygen Concentrator

Koncentrator tlenu

Table of Content

Spis treści

English	2-34
Polski	35-67

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Updated user manual

Using the QR code, download the current user manual containing the latest, updated information on the use, maintenance and troubleshooting of the device. The manual is available at welltec.pro/r/manual/ocm

Dear Customer,

Thank you for your trust and for choosing the **Welltec** Oxygen Concentrator. We hope that it will provide you a clean oxygen for many years.

Table of Content

Contents of the package	3
Product Introduction	4
First use	5
Precautions and warnings.....	6-9
Safety tips for oxygen Inhale	10
Safety tips for Atomization.....	11
Symbols	12
Usage	13-15
Remote controller	16
Atomization operation.....	17
Usage of Headset	18
Care and Cleaning.....	19-21
Specification	22
Troubleshooting	23-26
Pulse Oximeter introduction.....	27
Pulse Oximeter precautions and warnings	27-28
Pulse Oximeter information.....	29
Usage of Pulse Oximeter.....	30
Pulse Oximeter Battery Charging	31
Pulse Oximeter care and cleaning.....	31-32
Pulse Oximeter specification	33
Pulse Oximeter troubleshooting.....	34

Contents of the package

Oxygen Concentrator

remote control

pulse oximeter

2x nasal oxygen tube 1,6m

2x nasal oxygen tube 6m

2x oxygen mask

2x nebulizer kits

2x oxygen headset

splitter for connecting 2 people

4x sponge filter

2x air filter

pulse oximeter usb cable

Product Introduction

First use

NOTE. Before start check the concentrator has been connected with power supply and make sure it in power off status.

1. Remove the top cover part of humidification bottle and fill pure water (or distilled water), make sure the water level is between the „MAX“ and „MIN“ mark of the bottle and then cover the top part.

2. Install the humidification bottle on the humidifier bottle holder (A) of concentrator and secure it with band (B).

3. Connect the humidification bottle (A) to the oxygen outlet (C) using a connecting tube (B) and connect the other port of the humidification bottle (D) to the nasal oxygen tube (as shown below).

4. Press the power switch to position „I“ (ON). The power is on when the power indicator turns green. The oxygen concentrator is in working state.

5. Wear a nasal oxygen tube, or oxygen mask or oxygen headset.

6. Press the button  on the front panel to start oxygen supply.

7. Adjust the cannula to appropriate position for easy inhale of the pure oxygen and get maximum comfort. The oxygen density could be reached to 90% within **12 minutes**.

8. Turn off the equipment after use. Pull out the power pin and take care of the cannulas and oxygen tubes cleanliness for another use.

9. If no use of concentrator for a long-term, the water in humidifier bottle should be drained. Save the humidifier bottle after cleaned and dried.

Warning!

Every few seconds the concentrator makes „TAPOO“ sound while operating, it is the normal ventilation sound.

Oxygen uptake time and oxygen flow adjustment should follow the physician's advice.

Check that the tubing or accessory is clogged, kinked or the humidifier is damaged if the flow value is lower than 0.5L/min.

Precautions and warnings

To prevent injury or damage to property, read the following instructions.

Risk of health loss or death

Risk of health loss. Hazardous Substances

Prohibited

Should be done

To prevent oxygen concentrator's possible failures or power shut down, people in urgent need of oxygen and seriously ill patients must prepare other oxygen-supply devices for emergency use (etc: oxygen cylinders, oxygen bags).

The concentrator is suitable for oxygen supplement, and it is not intended to be life supporting or life sustaining.

Do not use the equipment near flammable materials as grease oil, detergent ect. In a certain pressure, oil, grease or grease substances which contacted with oxygen will self-ignition and intense combustion. These substances must be kept away from oxygen concentrator, piping, connectors, and all other oxygen devices. Do not use any lubricant other than the manufacturer's recommendation.

Oxygen is combustion-supporting gas. No smoking near the working oxygen concentrator. The oxygen concentrator should be kept away from matches, burning cigarettes and other objects of high temperature or fire. Textiles and other normally non-combustible materials are easily ignited and boiled in oxygen-enriched air. Ignorance of this warning may result in serious fire, property damage, and personal injury or death.

Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions, if the oxygen concentrator is turned on, but not in use; the oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.

Precautions and warnings

If any object or liquid enters the unit, disconnect the power plugs immediately, and have them tested by the professional person before re-use.

Do not disassemble.

Check the electrical connection regularly.

Use only the original plug without an extension cord.

Make sure that the device is properly connected to the power outlet.

50cm or more should be left between concentrator and the walls, windows, furniture and other similar objects

Ensure there is no humidification devices in the same room or within 2 meters around while using the concentrator.

Do not allow the device to get wet.

Keep away from running water.

Do not use in an explosive and corrosive atmosphere.

Do not use near open flame.

Do not place oxygen concentrator in parallel or in series with other concentrators or oxygen treatment equipment.

Do not block the air inlet / outlet under no circumstances.

Turn off the device before disconnecting the power plug.

After cleaning filter components, they must be totally dried before re-use.

Do not expose device to direct sunlight, wind or rain.

Do not place next to a heat source or steam.

Before turning on, check that the supply voltage from the rating plate matches the voltage of the local grid.

Do not pierce or set fire to the device

Do not spray the device with insecticide.

Do not touch with wet hands.

Do not clean with alcohol or solvents.

Precautions and warnings

Be careful to transport or convey, don't let it reversal or thwart, tilt angle not greater than 5°.

Geriatric, paediatrics or any other patient unable to communicate discomfort can require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to the responsible care giver to avoid harm.

Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.

Do not touch the cabinet or air outlet with your hands during the operation of the oxygen concentrator to avoid injury caused by overheat.

The temperature near the air outlet will be higher while the concentrator is operating, do not close to the air outlet to avoid being burned.

Do not move the oxygen concentrator when charged.

Personal and family use should be in accordance with the guidance of a physician when the concentration is higher than 93%.

Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.

Do not lubricate fittings, connections, nasal cannula, or other accessories of the oxygen concentrator to avoid the risk of fire and burns.

Open flames during oxygen therapy are dangerous and is likely to result in fire or death. Do not allow open flames within 2 m of the oxygen concentrator or any oxygen carrying accessories.

The device should be placed on a flat and stable surface to prevent water leakage.

Do not tread, sit on or lie on the oxygen concentrator.

Precautions and warnings

Wait 5 minutes before restarting the device.

Always grasp the plug when disconnecting the device.

Call the emergency hot line and seek the help of professional health care workers immediately if any discomfort is felt or accident happened while using the concentrator.

Attention!

If any of the following events occur, turn off the device immediately:

- One or more switches do not work.
- A short circuit has occurred.
- Disturbing overheating of the power cord / plug.
- Burning smell or disturbing sounds and vibrations.
- Other malfunctions.

If the power cord is damaged or destroyed, it must be replaced by the manufacturer, service or another authorized entity.

The device can be used by people with limited physical, sensory and mental performance on condition of assurance appropriate supervision and instruction in the field of use. Children should not play or use the device.

Attention!

The device is intended for use in rooms with a temperature from 5 to 40 ° C, humidity below 75%, barometric pressure range from 860 hPa to 1060 hPa and altitude below 4000 m.a.s.l. Operation of the device outside the specified temperatures, humidity, altitude, and atmospheric pressures may reduce the oxygen concentration.

Attention!

When using the concentrator, care shall be taken to prevent strangulation due to cables and hoses, particularly due to excessive length.

Safety tips for oxygen Inhale

No smoking while using the oxygen concentrator.

Follow the physician's guidance for it is used for medical treatment.

Do not move the oxygen concentrator when connected.

Patients with plastics allergy or vulnerable to oxygen poisoning are prohibited from using the device.

Clean and replace the air filter in case of the block of the outlet and outlet of oxygen delivery and affect the life of oxygen concentrator.

In order to keep the concentrator working normally, please substitute the filter sponge both in outlet and delivery with those provided or certified by manufacturer.

The oxygen flow rate shouldn't be too high, or upon the request of the physician. If the person or service person feels that the oxygen content is insufficient, please contact with doctor immediately and adjust the flow according to the doctor's instructions; patients with severe lung disease should consult professional physician for flow level.

The atomization device and the nasal cannula as accessories provided by our company is the test sample. The size of the atomization cup is 6-12ml, and the length of the nasal cannula is 1.6m and 6m. Use the nasal cannula personally to keep it clean and health; clean and disinfect it frequently; It is suggested using atomization device and disposal nasal cannula produced by qualified manufacturers.

The water in the bottle should not be too much in case of the overflow (please keep the water level between MAX and MIN), and change the water often. Different bottles may affect the performance of the oxygen concentrate, please replace the bottle with the original bottle provided by manufacturer.

Excessive oxygen inhalation can cause some damage, e.g. oxygen poisoning including CO₂ retention, newborns blind due to excessive oxygen intake, irritating dry cough, nausea, vomiting and headache, nasal duct injury, nasal bleeding, please use the oxygen concentrator with caution.

Safety tips for Atomization

 Check the connection of atomization assembly and oxygen concentrator is reliable before using. And check the leakage of atomization assembly, leakage will affect the atomization of the liquid.

 The water added in atomizing cup should not exceed 12ml each time.

 After the atomization assembly is connected, the flow level of the oxygen generator should be adjusted to the minimum.

The water in the atomizing cup must be clear, no precipitation and impurities, otherwise the straw and shunt will be blocked.

To prevent cross infection, it is recommended that each patient use different atomization assembly.


Different atomization assembly may affect the atomization effect, please use the original atomization assembly provided or certified atomization assembly by manufacturer.

Symbols

Some symbols may appear on your equipment.

Symbol	Description	Symbol	Description
	General warning sign		No open flame: Fire, open ignition source and smoking prohibited
	Caution		No smoking
	Keep dry		Fragile-handle with care
	This side up		Increase or decrease (Knob)
	ON (power)		Class II equipment
	OFF (power)		Alternating current
	Type BF applied part		Date of manufacturer
	Refer to instruction manual/booklet		Manufacturer
	Consult instructions for use		Serial number
	Degree of protection provided by enclosure		Batch code
	Part Number		Fuse
	Max. Stack Quantity		WEEE symbol
	CE mark		Authorized representative in the European community

1. Power ON/OFF



Press the I/O power switch on the front of concentrator to position „I“ (ON). The oxygen concentrator is in standby mode. Press the  button to turn on / off the device.

In the three seconds before oxygen concentrator restarts, the system will be automatically on self-check during which time the buzz and LCD without light belongs to normal situation. Oxygen concentration can reach a stable status after running about 12 minutes of concentrator.

2. 5. Timer

Press the  or  button to program the time after which the device will turn off. The possible settings are 1 - 10 hours. Then press button  to confirm the setting. When the shut hours are set, the system comes into count down time and the oxygen concentrator LCD (5) will show the remaining time. When remaining time becomes 0, the oxygen concentrator will go to standby status.

3. Fault code

The error code will be displayed in case of concentrator failure. Check the fault code table on page 26.

4. Total Running Time

The total running time is the sum of the working hours of the oxygen concentrator starting from the first use.

88888_H

5. Time / Timer

The Time / Timer display shows how long the concentrator works from its start-up or shows the remaining time until switching off in case of setting the Timer.

88:88

6. Power indicator

This indicator lights green when the oxygen concentrator is connected to the main power supply.

I/O

7. Air pressure / Concentration / Temperature / O2 sensor / Low voltage indicator

- This indicator lights yellow if the system air pressure is higher than 260kPa or lower than 20kPa and the screen shows the corresponding error code. Alert with buzzer sound.
- When the oxygen concentration is lower than 82%, the indicator lights up in yellow. Alert with buzzer sound.
- When the temperature of the system exceeds the maximum allowable temperature limit (70°C), the oxygen generator stops working. The indicator lights yellow and the screen shows the corresponding error code. Alert with buzzer sound.
- After tuning on, the oxygen concentrator does not receive any information from the O2 sensor. The indicator lights yellow and the screen show the corresponding error code. Alert with buzzer sound.
- When the power supply voltage is lower than 85% of the standard voltage, the indicator lights yellow and the screen shows the corresponding error code. And the oxygen concentrator alerts with buzzer sound.

A table with the individual error codes can be found on page 26.


8. Power failure alarm indicator

Accidental power off happened while operating the concentrator, this indicator lights red and flashes with buzzer alarm tone. This can prompt the user to disconnect the power and check the power supply system.


Flow knob

Rotate the knob on the front panel to adjust the flow. Rotate the knob to the left to increase the flow and to the right to decrease the flow. The greater the flow value, the greater the air flow, oxygen purity will be reduced. The recommended maximum flow rate of the OCM10 series oxygen concentrator is 10L/min. The recommended maximum concentration is $90\% \pm 3\%$.

Function between Oxygen concentration of concentrator and rated flow (outlet nominal pressure is 0kPa)

The effect of altitude on oxygen concentration delivered by oxygen concentrator:

Altitude(m)	≤800	800 - 1500	1500 - 3000	3000- 4000
OC0M10	≥90%	≥80%	≥72%	≥63%

Remote controller

Overview:

1. Working indicator (the indicator lights green when using buttons of remote controller)
2. Timing 1 hour button
3. Start/Stop oxygen button
4. Timing 2 hours button
5. Timing 3 hours button

Press button ☺ on the remote controller when the concentrator power is on and in standby status.

Press the button ☺ while the concentrator operating to stop the oxygen supply.

Press timer button on the remote controller to set timer length. When the remote controller cannot work, the battery may not have sufficient power. Please replace the battery periodically with the same battery and be careful of the affirmative and negative pole.

Warning!

- The concentrator should be in the remote range when using the controller.
- The maximum remote control distance is 5m.
- Please take the batteries out of the controller if unused for a long time.

Atomization operation

1. Press the power switch to position „I“ (ON) to turn on the concentrator.
2. Connect the atomized mask or mouthpiece to the medicine cup.
3. Unscrew the nut of the atomizer outlet; connect the atomization tubing to the outlet.
4. The press the button  on the front panel, and using knob to adjust the flow to minimum.
5. Put the atomization mouthpiece into the mouth or wear the mask and start the atomization therapy (follow the doctor's instructions for the atomization time).
6. After the atomization is finished, press the button  to stop atomizing. Remove tubing and turn off the oxygen concentrator by pressing the I/O power switch.


Atomization performance

The atomization outlet pressure is 60kPa to 250kPa, and the maximum air flow rate is $\geq 3\text{L/min}$. The atomization rate of the atomization assembly is $\geq 0.2\text{ mL/min}$.

Warning!

- Oxygen uptake and atomization cannot be used at the same time.
- Adjust the flow rate to minimum value when using the atomization function.

Usage and cleaning of Headset

Overview:

- | | |
|---------------------|-----------------------------|
| 1. Snap hook | 4. Fast hole |
| 2. $\triangle$ mark | 5. Top cover |
| 3. Slow hole | 6. Outlet transition sleeve |

Usage of headset

Headset can be worn on head, it makes user comfortable with its excellent personal feeling. Unique perfuming function: add different essential oil into liquid cup, this aromatherapy can be combined with oxygen at the same time.

Adjust the flow rate of headset

Headset has two levels flow rate: fast and slow. Turn the top cover (5) of Headset and make the „ $\triangle$ “ mark alignment (2). „Click“ sound will be heard while rotating the headset.

Wearing the Headset

Wear the Headset according to the figure shown during oxygen inhalation.

Usage of liquid cup

Add aromatic oil: turn the top cover (1) according to the arrow „OPEN  CLOSE“ (2) to open the cup, add the proper amount of aromatic oil, and then cover the top cover to the CLOSE position. Install the liquid cup correctly and then you can turn on the concentrator.

Care and Cleaning

Manufacturers recommend that the oxygen concentrator's running time is not less than 30 minutes each time. Do not turn on or off the oxygen concentrator frequently. Turn on the concentrator at least after 5 minutes after turning off to protect the compressor's life.

Warning!

Before cleaning the device, disconnect it from the power supply.

Care and Cleaning of cabinet

It's helpful for normal use to keep the oxygen concentration clean for long. The outside of the cabinet is cleaned at least once a month. Disconnect the equipment from the power supply before cleaning.

Concentrator cabinet cannot be washed with water, you can wipe its surface with a clean soft cotton or sponge. Do not pore the liquid into the cabinet gap. A mild home neutral detergent can be used to wipe the cabinet.

Care and Cleaning of humidifier bottle

Replace the water in bottle every day. Empty the humidifier each time after using the equipment.

Use clean water to clean the humidifier bottle. If there is any besmirch, you can use warm neutral detergent or solution of white vinegar and hot water with proportion of 1:10 to wash it. Clean the bottle with hot water, and then fill it with distilled water within the specified water level before using.


Care and Cleaning

Air filter

Cleaning and changing the air filter is important to protect and extend the life of the compressor and molecular sieve of the concentrator. It is recommended to clean the filter after about 100 hours of concentrator operation and. If the filter is damaged or worn out, replace it with a new one.

Disassembly air filter

The filter is located on the back of the concentrator. Remove the cover (3) by pressing the two buckles down and pulling outward, then remove the filter cover (2), and pull out the air filter (1).

Cleaning air filter

Clean the filter with clean water. If the filter is too dirty, you can use mild soap or detergent, then rinse thoroughly and air dry it. The filter must be completely dry before it can be installed on the oxygen concentrator.


Warning!

Do not operate oxygen concentrator without filter installation or with wet filter. Otherwise, it will damage the oxygen concentrator.

Care and Cleaning

Overload protection

When you suspect or determine the overload protection (turn on the main unit when the power connection is normal, a power failure alarm occurs), you can remove the cover (2) on the back of the concentrator and press the Reset (1) button.

Cleaning Headset

Rinse or immerse the disassembling liquid cup in warm water for about 15 minutes. Add appropriate amount of vinegar in the water for cleaning more hygienic. Clean the liquid cup with diluted cleaning agent, and then wash it by water.

Do not wash the above components with boiled water to prevent heat distortion. All components must be dried and stored after cleaning. Use clean soft cloth dipped in medical alcohol to wipe the Headset outlet.

Specification

Welltec Model	OCM10
Factory Model	KSOC-10
Rated flow	10l/min
Max oxygen concentration (V / V%)	90% ± 3%
Operating mode	Continuous
Control	• Control Panel • Remote control
Display	Yes
Functions	• Timer • Atomization • Total running time • Error display
Atomization	≥0.2 l/min
Atomizing port pressure range of concentrator with atomizing	60kPa - 250kPa
Atomization outlet pressure	30kPa -80kPa
Equipment	pulse oximeter 2x nasal oxygen tube 1,6m 2x nasal oxygen tube 6m 2x oxygen mask 2x nebulizer kits 2x oxygen headset splitter for connecting 2 people 4x sponge filter 2x air filter
Working range - temperature	from 5 to 40 °C
Working range - humidity	up to 75 %
Liquid protection level	IP21
Storage - temperature	from -40°C to +55°C
Storage - humidity	≤93%
Fuse	F8AL250V
Noise level	60 dB(A)
Power	820 VA
Power supply	220-230 V, 50 Hz
Weight	23 kg
Dimensions	380 x 320 x 620 mm
Package size	465 x 390 x 700 mm

Oxygen Concentrator Alarm and Indicator System

Alarm system design aims at monitoring the working oxygen concentrator in case of such situations as power off, abnormal pressure or indicator of running condition of the equipment. All alarms of this equipment are technical alarms.

Alarm	Alarm reason	Visual	Possible solutions
Power off, Triple + double + triple + double beep	The network power has been disconnected during operation.	Flashing frequency of the red alarm lamp: 1.4Hz - 2.8Hz.	Please turn off the power immediately. If an alarm still exists after confirming that the power supply and connection are normal, please turn off the oxygen concentrator and contact the place where you bought the device.
Low concentration, Single beep	Oxygen concentration is less than 82%.	Alarm indicator lights yellow.	Continue to use and contact the seller. Spare oxygen should be prepared for those who are in urgent need of oxygen.
Abnormal air pressure, Single beep	The internal tubing pressure of the concentrator is higher than 260kPa or lower than 20kPa.	Alarm indicator lights yellow and the screen displays error code „E05” or „E02”.	Please shut down the concentrator in the case of pressure alarm, and check and make sure the cleanness of outlet and inlet of air without any blocks and then restart the equipment. Please turn off the oxygen concentrator and inform the distributor if the alarm is still on.

If the problem cannot be solved, contact the place where you bought the device.

Troubleshooting

Alarm	Alarm reason	Visual	Possible solutions
High Temperature, Single beep	The oxygen concentrator system temperature is higher than the maximum allowable temperature limit.	The indicator lights yellow and the screen displays code "E35".	Turn off the power immediately. Check and make sure the cleanness of outlet and inlet of air without any blocks and then restart the Equipment. Please turn off the oxygen concentrator and inform the distributor if the alarm is still on.
Oxygen sensor communication failure, Single beep	The oxygen concentrator did not receive signal from the oxygen sensor.	The indicator lights yellow and the screen displays code "E31".	Contact the seller. Spare oxygen should be prepared for those who are in urgent need of oxygen.
Low Voltage, Single beep	The supply voltage of the concentrator is 85% lower than the standard voltage.	The indicator lights yellow and the screen displays code "E03".	Turn off the concentrator immediately, and restart it after confirming that the voltage of the network power supply is normal.

If the problem cannot be solved, contact the place where you bought the device.

Troubleshooting

Before reporting a bug, check the solutions below.

Problem	Possible Causes	Possible solutions
After turned on the switch, the light, alarm system and oxygen concentrator don't work.	• The plug is not properly connected to power outlet. • No power. • Fuse is broken.	• Insert the plug to the socket correctly. • Check the power supply. • Replace the fuse.
After turned on the switch, the light works, but the oxygen concentrator doesn't work.	• Air compressor protection. • Inlet or outlet jammed. • Ambient temperature is lower than 5°C	• If oxygen concentrator shutdown after restarting of 45 minutes, please contact the seller. • Clean the filter. Check if anything jammed intake. Check if objects blocking the cooling. • Rise the ambient temperature.
Atomization amount is insufficient or cannot be atomized while using the atomization.	• The atomizing assembly installed improperly or incorrectly. • The atomizing assembly outlet is blocked.	• Install the assembly properly and correctly. • Replace the atomizing assembly If this condition still exists, stop using the atomization assembly and contact the seller.
Cannot obtain the requested current capacity.	• Nasal oxygen tube is blocked or damaged. • Mask jammed or damaged • Humidifier bottle is blocked or damaged. • Oxygen tube surpasses the stipulation or has the bend.	• If the flow is normal, take down the nasal tube, clean, correct the dead bend or replace it. • If the flow is normal, take down the humidifier bottle, clean or replace it.

Troubleshooting

Before reporting a bug, check the solutions below.

Problem	Possible Causes	Possible solutions
Remote controller cannot work.	• The concentrator is out of the remote range. • Lower battery capacity or no power.	• Using remote controller within the remote range. • Replace battery.

Fault code table

Code	Fault Description
E02	During operation, the pressure suddenly drops below the limit pressure (20kPa).
E03	The supply voltage of the concentrator is lower than 85% of the standard voltage
E05	The pressure exceeds the limit pressure (260kPa) during operation.
E31	Cannot receive data from oxygen sensor.
E35	Compressor temperature control switch is not connected or exceeds the allowable range.

If the problem cannot be solved, contact the place where you bought the device. Non-professional maintenance personnel or personnel without our authorization are strictly prohibited to open the concentrator cabinet for maintenance.

Check System Gas Leakage and Gas Flowrate

Connect the nasal cannula to the gas outlet connector of the oxygen concentrator or, if used, to the bubble humidifier outlet connector per the manufacturer's instructions.

With the oxygen concentrator turned on adjust the flowmeter to the desired flowrate. Gas should be flowing freely to the nasal cannula. You should be able to hear or feel the flow of gas to the prongs of the nasal cannula.

Wave your hand in front of the uprongs. If you do not feel the gas flowing, check the cannula connections for leaks. or place the end of the nasal cannula under the surface of a half full cup of water and look for the bubbles.

Pulse Oximeter introduction

Overview:

1. Power On/Off switch:

- Turn on the oximeter: Press the ON/OFF button to turn on the device.
- Turn off the oximeter: Press this button for 2 seconds to turn the device off.
- Press this button to highlight the display when the display is in low brightness.

2. Bluetooth icon

- Visible: means that the oximeter is connected to a Bluetooth device.
- Not visible: means that no device is currently connected via Bluetooth.

3. Heart symbol

4. PR value: Pulse rate per minute

5. PR value: Pulse rate per minute

6. Saturation value (SpO2).

7. Battery capacity Indicator

-  - filling indicates battery level.
-  - the battery is discharged.

8. Micro USB connector: charging battery with specified charging cable connected to adaptor.

9. Rubber finger cover

Pulse Oximeter precautions and warnings

This device is not intended for medical treatment.

This device is not intended to use by infants and neonates.

Do not disassemble.

Do not use this device with the defibrillator or other high-frequency equipment.

Explosion hazard: Do not use this device in an explosive or corrosive atmosphere.

The measuring data displayed on the device are for reference only and cannot be directly used for diagnostic interpretation.

Pulse Oximeter precautions and warnings

Carefully read this manual about all safety information, operation and specifications before using this oximeter.

In some circumstances, the device may interpret motion as good pulse quality. Minimize patient motion as much as possible.

The presence of high ambient light may cause inaccurate SpO₂ measurements.

Do not connect it to any external device or operate it for measurement when this device is being charged.

Do not use this device in situations where alarms required. The oximeter does not support alarms.

Users who are allergic to rubber cannot use this product.

Check and change the applied site according to the different circumstances of the user while using this device for long-time continuous monitoring. It is recommended to check the finger after 2 hours. Change other finger to be measured if the finger seems to be unusual. If the oximeter is too tight because the application site is too large or becomes too large due to edema, excessive pressure for prolonged periods may result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.

Do not allow the device to get wet.

Do not pierce or set fire to the device.

Keep away from running water.

Do not spray the device with insecticide.

Do not use in an explosive and corrosive atmosphere.

Do not touch with wet hands.

Do not use near open flame.

Do not clean with alcohol or solvents.

Do not use this device near by the source of electromagnetic interference, such as mobile phones or radio transmitter.

Do not use this device if the device cannot achieve satisfactory results.

Avoid placing the sensor on extremities with an arterial catheter, or intravascular venous infusion line.

Always install or carry the device properly to avoid damage caused by drop, impact, strong vibration or other mechanical force.

Check the device before using. Do not use it if there is significant damage.

Do not use the functional tester to access the SpO2 accuracy.

Pulse Oximeter information

Pulse oximeter using spectrophotometry, by detecting the blood red and infrared absorption of light to obtain the oxygen saturation (SpO₂) and pulse rate (PR).

Oxygen saturation is a term referring to the fraction of oxygen-saturated hemoglobin (HbO₂) relative to total hemoglobin (Hb) in the blood, that is, the oxygen saturation in the blood, which is an important physiological parameter of the respiratory and circulatory system. Many respiratory diseases can cause oxygen saturation decrease. The body's automatic regulatory dysfunction caused by anesthesia, major surgery trauma, and injury caused by some medical examination etc., are likely to lead to the oxygen saturation decrease, resulting in patients with dizziness, weakness, vomiting and other symptoms; severe cases will be life-threatening. SpO₂ has important significance in the field of clinical medical. And timely understanding of the patient's blood oxygen saturation situation will help physician find problems.

Intended Use and Contraindication

The fingertip pulse oximeter is intended to measure the pulse oxygen saturation of arterial hemoglobin and pulse rate of adults and pediatrics in home care and medical outpatient environment.

Contraindications: users allergic to rubber cannot use this product.

Applications:

This device is suitable for the following people:

- People with vascular disease, such as: coronary heart disease, hypertension, hyperlipidemia, cerebral thrombosis;
- People with respiratory diseases, such as: asthma, trachitis, chronic bronchitis, chronic cor pulmonale, chronic obstructive pulmonary disease;
- Old people above 60 years old;
- People who works more than 12 hours a day;
- People works on extreme exercise or under the alpine hypoxic environment;
- Long-term alcohol people.

Usage of Pulse Oximeter

1. Insert finger to oximeter's cover fully to obtain good measurements.
2. Press the ON/OFF button to turn the oximeter on.
3. SpO2 and pulse rate value displays on screen.

Notes:

- **The oximeter cannot be measured while charging.**
- Finger cannot apply nail polish and other cosmetics.
- Do not shake the finger, hand or body during the measuring.
- Tape or other light obstructions around the applied site may affect the accuracy of SpO2 and pulse rate.
- The oximeter will shut down if no any actions within 10 seconds after it turning on.

The following factors may influence the accuracy of measurements:

- Ambient light;
- Physical movement (patient or imposed motion);
- Diagnostic testing;
- Low perfusion;
- Electromagnetic interference, such as MRI environment;
- Electrosurgical units;
- Dysfunctional haemoglobin, such as carboxyhemoglobin(COHb) and methemoglobin(MetHb);
- The oximeter placed improperly;
- Drop of arterial blood flow to immeasurable level caused by shock, anemia, low temperature or vasoconstrictor;
- Dopamine, procaine, propranolol, lidocaine, bupivacaine and other drugs can cause more severe bias of SpO2
- For anemia hypoxia and toxic hypoxia, the SpO2 value is only for reference, as some patients with severe anemia still show higher SpO2 value.

Pulse Oximeter Battery Charging

The oximeter uses a built-in non-removable rechargeable lithium battery. Symbol  will flash when battery capacity is low. Charge the battery according to the following steps.

1. Take out the cable and connect it to oximeter's Micro USB connector.
2. Plug the other end of cable to USB adaptor. Battery symbol appear and flashes on the screen while charging.

The performance of the rechargeable lithium battery may deteriorate over time. If the operating time of the battery is noticeably shorter than that stated in the specifications, contact the place where you bought the device.

Pulse Oximeter care and cleaning

Use only the substances approved by us and methods listed in this chapter to clean or disinfect the oximeter.

We make no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. For the method to control infection, consult your hospital's Infection Control Officer or Epidemiologist.

To avoid damage to the oximeter, follow these rules:

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse the oximeter into liquid.
- Do not pour liquid onto the oximeter or accessories.
- Do not allow liquid to enter the oximeter case.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).

Pulse Oximeter care and cleaning

Recommended cleaning agents:

- Ethanol (75%)
- Isopropanol (70%)

Warning!

Turn off the oximeter before cleaning and unplug the charging cable.

Do not disinfect the oximeter with high temperature, high pressure, gas vapor or liquid immersion. Clean and disinfect the oxygen meter according to the manufacturer's requirements.

All servicing should be done by qualified professionals; user or other person cannot disassembly or repair the oximeter.

If you spill liquid on the device or accessories, wipe it clean immediately. If the oximeter does not work properly, contact the place where you bought the device.

Cleaning

Follow the below steps to clean the oximeter:

- 1.** Make sure turning off the oximeter.
- 2.** Clean the exterior surface of the oximeter with a soft cloth dampened with a cleaning agent.
- 3.** Wipe off all the cleaning solution with a dry cloth after cleaning if necessary.
- 4.** Dry the oximeter in a ventilated, cool place.

Disinfection

Disinfection may cause damage to the oximeter. We recommend that you disinfect the oximeter only when necessary. Clean the oximeter before disinfection. Recommended cleaning agent is Ethanol (75%) The recommended times of disinfection is three hundred.

Note: Never use ETO or formaldehyde for disinfection.

Pulse Oximeter specification

Model		F01LT
Operating mode		Continuous
SpO2 Saturation (SpO2%)	Measuring range	0% ~ 100%
	Measurement precision	70% ~ 100%, $\pm 2\%$
	Measurement resolution	1%
Pulse rate (PR)	Measuring range	18bpm ~ 250bpm
	Measurement precision	25bpm ~ 250bpm, ± 3 bpm
	Measurement resolution	1bpm
	Refreshing rate	1s
Control		• Control Panel • WiFi
Display		Yes - LCD
Transducer wavelength		Red light: 660 ± 3 nm; Infrared light: 940 ± 5 nm
Transducer maximum photic output consumption:		<30mW
Working range - temperature		from 5 to 40 °C
Working range - humidity		from 30% to 80%
Recommended max charge and discharge ambient temperature		Charge: 35°C, discharge: 45°C
Liquid protection level		IPX1, non AP/APG type
Storage - temperature		from -20°C to +60°C
Storage - humidity		from 10% to 95%
Battery capacity		300mAh
Charging voltage and current		5V ~ 160mA
Charging connector		micro USB connector
Charge time		2.5 hours (charge to 90% of capacity)
Run time		20 hours (Fully-charged new battery, ambient temperature: 25°C)
Shutdown delay		max 5 minutes after the low battery prompt first occurs
Weight		< 50g
Dimensions		72 x 32 x 27 mm

Note: The information about wavelength range and maximum photic output consumption can be especially useful to clinicians, for example, clinicians performing photodynamic therapy.


Pulse Oximeter troubleshooting

Before reporting a bug, check the solutions below.


Problem	Possible Causes	Possible solutions
The oximeter cannot start.	• Lower battery capacity or no power, • Oximeter might be damaged.	• Charge the battery. • Check the power supply.
SpO ₂ or PR value cannot be shown normally.	• Finger might not be inserted properly, • High levels of ambient light, • User is in low perfusion or user's oxyhemoglobin is too low to be measured, • Oximeter might be damaged.	• Do not use the oximeter in the environment with high ambient light, • Go to a hospital for diagnosis.
SpO ₂ or PR value is unstable.	• Finger might not be inserted properly, • Finger is trembling or patient is moving.	• Insert finger properly. • Please keep still.
The display is suddenly off.	• No finger inserted after turning on the oximeter, or it automatically shut down caused by lowest battery capacity.	• Press ON/OFF button to turn it on or charge the battery.

Manufactured for:

Manufactured by:

Warning! Do not dispose of the device in the household waste!

In accordance with the guidelines of the European Directive 2012/19 / EU on old, used electrical and electronic devices and its application in national law, unusable electrical equipment has been eliminated and electronic must be collected separately and recycled in an environmentally sound manner.

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Jiangsu Konsung Bio-Medical Science And Technology Co.,Ltd.
Shanghai International Holding Corp. GmbH (Europe)